

Axin-2

Catalog No: #AP80880



Package Size: #AP80880-1 50ug #AP80880-2 100ug #AP80880-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Axin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Axin-like protein,Axil,Axis inhibition protein 2,Conductin
Accession No.	Uniprot ID:O88566
Uniprot	O88566
GeneID	12006;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Axin2

Note: This product is for in vitro research use only