

Basal cell adhesion molecule

Catalog No: #AP80883



Package Size: #AP80883-1 50ug #AP80883-2 100ug #AP80883-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Basal cell adhesion molecule
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	B-CAM cell surface glycoprotein,Lutheran antigen,CD239
Accession No.	Uniprot ID:Q9MZ08
Uniprot	Q9MZ08
GeneID	282862;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BCAM

Note: This product is for in vitro research use only