

Plasma membrane calcium-transporting ATPase 1

Catalog No: #AP80887



Package Size: #AP80887-1 50ug #AP80887-2 100ug #AP80887-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Plasma membrane calcium-transporting ATPase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PMCA1, Plasma membrane calcium ATPase isoform 1, Plasma membrane calcium pump isoform 1
Accession No.	Uniprot ID:P11505
Uniprot	P11505
GeneID	29598;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Atp2b1

Note: This product is for in vitro research use only