

Autophagy-related protein 9A

Catalog No: #AP80889



Package Size: #AP80889-1 50ug #AP80889-2 100ug #AP80889-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Autophagy-related protein 9A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	APG9-like 1,Autophagy protein 9,ATG9
Accession No.	Uniprot ID:Q3T904
Uniprot	Q3T904
GeneID	540482;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATG9A

Note: This product is for in vitro research use only