

## Myc box-dependent-interacting protein 1

Catalog No: #AP80895



Package Size: #AP80895-1 50ug #AP80895-2 100ug #AP80895-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myc box-dependent-interacting protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Amphiphysin II,Amphiphysin-like protein,Bridging integrator 1,Amph2,Amphl
Accession No.	Uniprot ID:O08839
Uniprot	O08839
GeneID	117028;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Bin1

Note: This product is for in vitro research use only