

B-cell lymphoma 3 protein

Catalog No: #AP80896



Package Size: #AP80896-1 50ug #AP80896-2 100ug #AP80896-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B-cell lymphoma 3 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BCL-3,Proto-oncogene BCL3,BCL4,D19S37
Accession No.	Uniprot ID:P20749
Uniprot	P20749
GeneID	602;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BCL3

Note: This product is for in vitro research use only