

Bone morphogenetic protein 3B

Catalog No: #AP80899



Package Size: #AP80899-1 50ug #AP80899-2 100ug #AP80899-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bone morphogenetic protein 3B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	BMP-3B,Bone-inducing protein,BIP,Growth/differentiation factor 10,GDF-10,Bmp3b
Accession No.	Uniprot ID:P55108
Uniprot	P55108
GeneID	79216;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gdf10

Note: This product is for in vitro research use only