

BCL2/adenovirus E1B 19 kDa protein-interacting protein 3

Catalog No: #AP80900

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80900-1 50ug #AP80900-2 100ug #AP80900-3 1mg

Description

Product Name	BCL2/adenovirus E1B 19 kDa protein-interacting protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NIP3
Accession No.	Uniprot ID:Q12983
Uniprot	Q12983
GeneID	664;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BNIP3

Note: This product is for in vitro research use only