

# A disintegrin and metalloproteinase with thrombospondin motifs 1

Catalog No: #AP80901

Package Size: #AP80901-1 50ug #AP80901-2 100ug #AP80901-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADAM-TS 1
Accession No.	Uniprot ID:P97857
Uniprot	P97857
GeneID	11504;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Adamts1

Note: This product is for in vitro research use only