

Activin receptor type-2A

Catalog No: #AP80902



Package Size: #AP80902-1 50ug #AP80902-2 100ug #AP80902-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Activin receptor type-2A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Activin receptor type IIA,ACTR-IIA,Actrii,Acvr2
Accession No.	Uniprot ID:P38444
Uniprot	P38444
GeneID	29263;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Acvr2a

Note: This product is for in vitro research use only