

Basal cell adhesion molecule

Catalog No: #AP80903



Package Size: #AP80903-1 50ug #AP80903-2 100ug #AP80903-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Basal cell adhesion molecule
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Auberger B antigen,B-CAM cell surface glycoprotein,F8/G253 antigen,Lutheran antigen,Lutheran blood group glycoprotein,CD239,LU,MSK19
Accession No.	Uniprot ID:P50895
Uniprot	P50895
GeneID	4059;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BCAM

Note: This product is for in vitro research use only