

S-arrestin

Catalog No: #AP80905



Package Size: #AP80905-1 50ug #AP80905-2 100ug #AP80905-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	S-arrestin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	48 kDa protein, Retinal S-antigen, S-AG, Rod photoreceptor arrestin
Accession No.	Uniprot ID:P08168
Uniprot	P08168
GeneID	280922;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SAG

Note: This product is for in vitro research use only