

## Sarcoplasmic/endoplasmic reticulum calcium ATPase

3



Catalog No: #AP80907

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP80907-1 50ug #AP80907-2 100ug #AP80907-3 1mg

## Description

Product Name	Sarcoplasmic/endoplasmic reticulum calcium ATPase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	ChkSERCA3, Calcium pump 3
Accession No.	Uniprot ID:Q9YGL9
Uniprot	Q9YGL9
GeneID	395707;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: ATP2A3

Note: This product is for in vitro research use only