

Cystatin C Mouse Monoclonal Antibody(7F11)

Catalog No: #44153



Package Size: #44153-1 50ul #44153-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Cystatin C Mouse Monoclonal Antibody(7F11)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7F11
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu
Specificity	Cystatin C Mouse Monoclonal Antibody(7F11) detects endogenous levels of CST3
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Cystatin C
Other Names	CST3
Accession No.	Swiss-Prot:P01034Gene ID:1471
Uniprot	P01034
GeneID	1471;
Calculated MW	14kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1:1000-1:2000
Immunohistochemistry: 1:100-1:200


Images



Western Blot analysis of Cystatin C protein using antibody diluted at 1:1000



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using Cystatin C Mouse mAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using Cystatin C Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only