

Carbonic anhydrase 1

Catalog No: #AP80910



Package Size: #AP80910-1 50ug #AP80910-2 100ug #AP80910-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carbonic anhydrase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Carbonate dehydratase I,Carbonic anhydrase B,CAB,Carbonic anhydrase I,CA-I
Accession No.	Uniprot ID:P00915
Uniprot	P00915
GeneID	759;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CA1

Note: This product is for in vitro research use only