

Alpha-S2-casein

Catalog No: #AP80912



Package Size: #AP80912-1 50ug #AP80912-2 100ug #AP80912-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-S2-casein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Accession No.	Uniprot ID:P39036
Uniprot	P39036
GeneID	445515;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CSN1S2

Note: This product is for in vitro research use only