

Cathepsin G

Catalog No: #AP80913



Package Size: #AP80913-1 50ug #AP80913-2 100ug #AP80913-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Cathepsin G |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Vimentin-specific protease,VSP |
| Accession No. | Uniprot ID:P28293 |
| Uniprot | P28293 |
| GeneID | 13035; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ctsg

Note: This product is for in vitro research use only