

A disintegrin and metalloproteinase with thrombospondin motifs 9

Catalog No: #AP80916

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80916-1 50ug #AP80916-2 100ug #AP80916-3 1mg

Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADAM-TS 9,KIAA1312
Accession No.	Uniprot ID:Q9P2N4
Uniprot	Q9P2N4
GeneID	56999;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADAMTS9

Note: This product is for in vitro research use only