

Copper-transporting ATPase 1

Catalog No: #AP80918



Package Size: #AP80918-1 50ug #AP80918-2 100ug #AP80918-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Copper-transporting ATPase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Copper pump 1, Menkes disease-associated protein homolog, Mnk
Accession No.	Uniprot ID:Q64430
Uniprot	Q64430
GeneID	11977;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Atp7a

Note: This product is for in vitro research use only