

5-azacytidine-induced protein 2

Catalog No: #AP80919



Package Size: #AP80919-1 50ug #AP80919-2 100ug #AP80919-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	5-azacytidine-induced protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NF-kappa-B-activating kinase-associated protein 1, Nak-associated protein 1, Az2, Nap1, Tbkp2
Accession No.	Uniprot ID:Q9QYP6
Uniprot	Q9QYP6
GeneID	27215;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Azi2

Note: This product is for in vitro research use only