

Sphingomyelin phosphodiesterase

Catalog No: #AP80922



Package Size: #AP80922-1 50ug #AP80922-2 100ug #AP80922-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--------------------------------------------------------|
| Product Name | Sphingomyelin phosphodiesterase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Acid sphingomyelinase,ASMase,Asm |
| Accession No. | Uniprot ID:Q04519 |
| Uniprot | Q04519 |
| GeneID | 20597; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Smpd1

Note: This product is for in vitro research use only