

# Sodium/potassium-transporting ATPase subunit alpha-1

Catalog No: #AP80930

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP80930-1 50ug #AP80930-2 100ug #AP80930-3 1mg

## Description

Product Name	Sodium/potassium-transporting ATPase subunit alpha-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Na(+)/K(+) ATPase alpha-1 subunit,Sodium pump subunit alpha-1
Accession No.	Uniprot ID:P50997
Uniprot	P50997
GeneID	403992;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ATP1A1

Note: This product is for in vitro research use only