

Autophagy protein 5

Catalog No: #AP80931



Package Size: #AP80931-1 50ug #AP80931-2 100ug #AP80931-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Autophagy protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	APG5-like, Apoptosis-specific protein, APG5L, ASP
Accession No.	Uniprot ID:Q9H1Y0
Uniprot	Q9H1Y0
GeneID	9474;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ATG5

Note: This product is for in vitro research use only