

Cyclic AMP-dependent transcription factor ATF-1

Catalog No: #AP80933



Package Size: #AP80933-1 50ug #AP80933-2 100ug #AP80933-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclic AMP-dependent transcription factor ATF-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	cAMP-dependent transcription factor ATF-1,Activating transcription factor 1,Protein TREB36
Accession No.	Uniprot ID:P18846
Uniprot	P18846
GeneID	466;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATF1

Note: This product is for in vitro research use only