

A disintegrin and metalloproteinase with thrombospondin motifs 4

Catalog No: #AP80934

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80934-1 50ug #AP80934-2 100ug #AP80934-3 1mg

Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ADAM-TS 4,ADMP-1,Aggrecanase-1
Accession No.	Uniprot ID:Q9TT93
Uniprot	Q9TT93
GeneID	286806;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADAMTS4

Note: This product is for in vitro research use only