

Sodium/potassium-transporting ATPase subunit alpha-1

Catalog No: #AP80935

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP80935-1 50ug #AP80935-2 100ug #AP80935-3 1mg

Description

Product Name	Sodium/potassium-transporting ATPase subunit alpha-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Na(+)/K(+) ATPase alpha-1 subunit,Sodium pump subunit alpha-1
Accession No.	Uniprot ID:Q8VDN2
Uniprot	Q8VDN2
GeneID	11928;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Atp1a1

Note: This product is for in vitro research use only