

ADP-ribosylation factor-like protein 13B

Catalog No: #AP80938



Package Size: #AP80938-1 50ug #AP80938-2 100ug #AP80938-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | ADP-ribosylation factor-like protein 13B |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | ADP-ribosylation factor-like protein 2-like 1,ARL2-like protein 1,Protein hennin,Arl2I1 |
| Accession No. | Uniprot ID:Q640N2 |
| Uniprot | Q640N2 |
| GeneID | 68146; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Arl13b

Note: This product is for in vitro research use only