

## Cystatin C Mouse Monoclonal Antibody(5A2)

Catalog No: #44156

Package Size: #44156-1 50ul #44156-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	Cystatin C Mouse Monoclonal Antibody(5A2)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	5A2
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu
Specificity	Cystatin C Mouse Monoclonal Antibody(5A2) detects endogenous levels of CST3
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Cystatin C
Other Names	CST3
Accession No.	Swiss-Prot:P01034Gene ID:1471
Uniprot	P01034
GeneID	1471;
Calculated MW	14kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1:1000-1:2000

Immunohistochemistry: 1:100-1:200

## Images



Western Blot analysis of Cystatin C protein using antibody diluted at 1:1000



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using Cystatin C Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only