

Beta-secretase 1

Catalog No: #AP80941



Package Size: #AP80941-1 50ug #AP80941-2 100ug #AP80941-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-secretase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Aspartyl protease 2,ASP2,Beta-site amyloid precursor protein cleaving enzyme 1,Beta-site APP cleaving enzyme 1,Memapsin-2,Membrane-associated aspartic protease 2,Bace
Accession No.	Uniprot ID:P56818
Uniprot	P56818
GeneID	23821;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bace1

Note: This product is for in vitro research use only