

Cyclic AMP-dependent transcription factor ATF-2

Catalog No: #AP80944



Package Size: #AP80944-1 50ug #AP80944-2 100ug #AP80944-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclic AMP-dependent transcription factor ATF-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	cAMP-dependent transcription factor ATF-2,Activating transcription factor 2
Accession No.	Uniprot ID:O93602
Uniprot	O93602
GeneID	395727;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATF2

Note: This product is for in vitro research use only