

Isoaspartyl peptidase/L-asparaginase

Catalog No: #AP80945



Package Size: #AP80945-1 50ug #AP80945-2 100ug #AP80945-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Isoaspartyl peptidase/L-asparaginase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Asparaginase-like protein 1,Asparaginase-like sperm autoantigen,Beta-aspartyl-peptidase,Glial asparaginase,Isoaspartyl dipeptidase,L-asparagine amidohydrolase,Alp,Gliap,Hiob
Accession No.	Uniprot ID:Q8VI04
Uniprot	Q8VI04
GeneID	246307;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Asrgl1

Note: This product is for in vitro research use only