

Myc box-dependent-interacting protein 1

Catalog No: #AP80946



Package Size: #AP80946-1 50ug #AP80946-2 100ug #AP80946-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myc box-dependent-interacting protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Amphiphysin II,Amphiphysin-like protein,Box-dependent myc-interacting protein 1,Bridging integrator 1,AMPHL
Accession No.	Uniprot ID:O00499
Uniprot	O00499
GeneID	274;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BIN1

Note: This product is for in vitro research use only