

B-cell linker protein

Catalog No: #AP80947



Package Size: #AP80947-1 50ug #AP80947-2 100ug #AP80947-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	B-cell linker protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	B-cell adapter SH2 domain-containing protein,Cytoplasmic adapter protein,Src homology 2 domain-containing leukocyte protein of 65 kDa,SLP-65,BASH
Accession No.	Uniprot ID:Q9YGC1
Uniprot	Q9YGC1
GeneID	395733;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BLNK

Note: This product is for in vitro research use only