

Atrial natriuretic peptide receptor 1

Catalog No: #AP80950



Package Size: #AP80950-1 50ug #AP80950-2 100ug #AP80950-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Atrial natriuretic peptide receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Atrial natriuretic peptide receptor type A, ANP-A, Guanylate cyclase A, GC-A, Npra
Accession No.	Uniprot ID:P18293
Uniprot	P18293
GeneID	18160;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Npr1

Note: This product is for in vitro research use only