

Beta-2-microglobulin

Catalog No: #AP80957



Package Size: #AP80957-1 50ug #AP80957-2 100ug #AP80957-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Beta-2-microglobulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Lactollin
Accession No.	Uniprot ID:P01888
Uniprot	P01888
GeneID	280729;783680;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: B2M

Note: This product is for in vitro research use only