

Cystatin C Mouse Monoclonal Antibody(3B12)

Catalog No: #44158

Package Size: #44158-1 50ul #44158-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | Cystatin C Mouse Monoclonal Antibody(3B12) |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | 3B12 |
| Isotype | IgG1 |
| Applications | WB |
| Species Reactivity | Hu |
| Specificity | Cystatin C Mouse Monoclonal Antibody(3B12) detects endogenous levels of CST3 |
| Immunogen Type | Recombinant protein |
| Immunogen Description | Recombinant Protein |
| Target Name | Cystatin C |
| Other Names | CST3 |
| Accession No. | Swiss-Prot:P01034Gene ID:1471 |
| Uniprot | P01034 |
| GeneID | 1471; |
| Calculated MW | 14kd |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

Western blotting: 1:1000-1:2000

Images



Western Blot analysis of Cystatin C protein using antibody diluted at 1:1000

Note: This product is for in vitro research use only