

Bone morphogenetic protein 1

Catalog No: #AP80961



Package Size: #AP80961-1 50ug #AP80961-2 100ug #AP80961-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bone morphogenetic protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BMP-1,Mammalian tolloid protein,mTld,Procollagen C-proteinase,PCP
Accession No.	Uniprot ID:P98063
Uniprot	P98063
GeneID	12153;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bmp1

Note: This product is for in vitro research use only