

Complement component 1 Q subcomponent-binding protein, mitochondrial

Catalog No: #AP80963

Orders: order@signalwayantibody.com

Package Size: #AP80963-1 50ug #AP80963-2 100ug #AP80963-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Complement component 1 Q subcomponent-binding protein, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q3T0B6
Uniprot	Q3T0B6
GeneID	518321;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: C1QBP

Note: This product is for in vitro research use only