

Transforming growth factor-beta-induced protein ig-h3

Catalog No: #AP80964



Package Size: #AP80964-1 50ug #AP80964-2 100ug #AP80964-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transforming growth factor-beta-induced protein ig-h3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Beta ig-h3,MP70,MP78,BIGH3
Accession No.	Uniprot ID:P55906
Uniprot	P55906
GeneID	539596;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TGFBI

Note: This product is for in vitro research use only