

Basigin

Catalog No: #AP80969



Package Size: #AP80969-1 50ug #AP80969-2 100ug #AP80969-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Basigin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Glycoprotein CE9,OX-47 antigen,CD147
Accession No.	Uniprot ID:P26453
Uniprot	P26453
GeneID	25246;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bsg

Note: This product is for in vitro research use only