

BH3-interacting domain death agonist

Catalog No: #AP80970



Package Size: #AP80970-1 50ug #AP80970-2 100ug #AP80970-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BH3-interacting domain death agonist
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	p22 BID,BID
Accession No.	Uniprot ID:P55957
Uniprot	P55957
GeneID	637;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BID

Note: This product is for in vitro research use only