

# Autophagy-related protein 9A

Catalog No: #AP80971



Package Size: #AP80971-1 50ug #AP80971-2 100ug #AP80971-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Autophagy-related protein 9A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	APG9-like 1,mATG9,APG9L1
Accession No.	Uniprot ID:Q7Z3C6
Uniprot	Q7Z3C6
GeneID	79065;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ATG9A

Note: This product is for in vitro research use only