

Plasma membrane calcium-transporting ATPase 2

Catalog No: #AP80973



Package Size: #AP80973-1 50ug #AP80973-2 100ug #AP80973-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Plasma membrane calcium-transporting ATPase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PMCA2, Plasma membrane calcium ATPase isoform 2, Plasma membrane calcium pump isoform 2, Pmca2
Accession No.	Uniprot ID:Q9R0K7
Uniprot	Q9R0K7
GeneID	11941;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Atp2b2

Note: This product is for in vitro research use only