

ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 2

Catalog No: #AP80975



Package Size: #AP80975-1 50ug #AP80975-2 100ug #AP80975-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADP-ribosyl cyclase 2, Antigen BP3, BP-3 alloantigen, Bone marrow stromal antigen 1, BST-1, Cyclic ADP-ribose hydrolase 2, cADPr hydrolase 2, Leukocyte antigen 65, Ly-65, CD157, Bp-3, Bp3, Ly65
Accession No.	Uniprot ID:Q64277
Uniprot	Q64277
GeneID	12182;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Bst1

Note: This product is for in vitro research use only