

B- and T-lymphocyte attenuator

Catalog No: #AP80976



Package Size: #AP80976-1 50ug #AP80976-2 100ug #AP80976-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B- and T-lymphocyte attenuator
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	B- and T-lymphocyte-associated protein,CD272
Accession No.	Uniprot ID:Q7TSA3
Uniprot	Q7TSA3
GeneID	208154;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Btla

Note: This product is for in vitro research use only