

Scavenger receptor cysteine-rich type 1 protein M130

Catalog No: #AP80977



Package Size: #AP80977-1 50ug #AP80977-2 100ug #AP80977-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Scavenger receptor cysteine-rich type 1 protein M130
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CD163,M130
Accession No.	Uniprot ID:Q2VLH6
Uniprot	Q2VLH6
GeneID	93671;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cd163

Note: This product is for in vitro research use only