

BPI fold-containing family A member 2

Catalog No: #AP80978



Package Size: #AP80978-1 50ug #AP80978-2 100ug #AP80978-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BPI fold-containing family A member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Parotid secretory protein,PSP,Short palate, lung and nasal epithelium carcinoma-associated protein 2,C20orf70,SPLUNC2,UNQ510/PRO1025
Accession No.	Uniprot ID:Q96DR5
Uniprot	Q96DR5
GeneID	140683;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BPIFA2

Note: This product is for in vitro research use only