

A disintegrin and metalloproteinase with thrombospondin motifs 5

Catalog No: #AP80979

Orders: order@signalwayantibody.com

Package Size: #AP80979-1 50ug #AP80979-2 100ug #AP80979-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADAM-TS 5,ADMP-2,Aggrecanase-2,Implantin
Accession No.	Uniprot ID:Q9R001
Uniprot	Q9R001
GeneID	23794;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adamts5

Note: This product is for in vitro research use only