

Beta-1,4-galactosyltransferase 5

Catalog No: #AP80981



Package Size: #AP80981-1 50ug #AP80981-2 100ug #AP80981-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-1,4-galactosyltransferase 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Beta-1,4-GalTase 5,Beta-1,4-GalT II,UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 5,UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 5,Bgt-5
Accession No.	Uniprot ID:Q9JMK0
Uniprot	Q9JMK0
GeneID	56336;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:B4galt5

Note: This product is for in vitro research use only