

# Sodium/potassium-transporting ATPase subunit alpha-3

Catalog No: #AP80982

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP80982-1 50ug #AP80982-2 100ug #AP80982-3 1mg

## Description

Product Name	Sodium/potassium-transporting ATPase subunit alpha-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Na(+)/K(+) ATPase alpha-3 subunit,Na(+)/K(+) ATPase alpha(III) subunit,Sodium pump subunit alpha-3
Accession No.	Uniprot ID:Q6PIC6
Uniprot	Q6PIC6
GeneID	232975;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Atp1a3

Note: This product is for in vitro research use only